# boso medicus family 4

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**(D)** 2

**User  
Instructions**

**(GB)** 38

**Mode  
d'emploi**

**(F)** 74

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**(I)** 110

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**(E)** 146

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# Table of contents

---

Product features.....	40
Contents of package .....	41
Quick guide .....	42
Introduction.....	43
Purpose.....	44
Blood pressure values.....	45
Systolic and diastolic blood pressure value .....	45
WHO blood pressure cut-off values .....	45
Starting up.....	47
Inserting batteries .....	47
Battery status.....	48
Set Time and Date.....	48
General instructions for self-measurement.....	50
Preparing to measure your blood pressure.....	52
Note on cuff sizes.....	52
Attaching the cuff.....	53

## Table of contents

---

Select an user.....	56
Measuring your blood pressure.....	56
Interrupting a measurement.....	57
Measurement display.....	58
Note on guest mode.....	58
Completing the blood pressure record card / Retrieving the average value and WHO classification.....	60
Displaying individual stored values.....	62
Mains operation.....	63
Error messages.....	64
Cleaning and Disinfection.....	66
Warranty Conditions.....	67
Customer service / Disposal.....	68
Technical data.....	69
Calibration checks - Testing instructions (for trained specialist staff only).....	72

# Product features

Battery status display

see page 48

Measurement in progress

see page 57

Irregular heartbeat during measurement

see page 58, 65

Movement during measurement

see page 57, 59

Rate of recurrence of cardiac arrhythmias

see page 60

Stored value

see page 61, 62

Select an user

see page 56

# Contents of package

1 **boso-medicus family 4**  
blood pressure monitor


4 LR 6 (AA) batteries

1 CA 04 universal cuff

1 User  
instructions

1 Case

1 Blood pressure  
record card

## Quick guide

- ➔ Insert batteries, ensure correct positioning (P. 47)
- ➔ Connect cuff to monitor (P. 52)
- ➔ Attach cuff (P. 53)
- ➔ Start (and stop) a measurement by pressing the start button (P. 56)

Measured systolic, diastolic and pulse values are displayed after the measurement is completed (P. 58)

- ➔ Retrieving values from memory (P. 60): Press the M-button with the monitor switched off. The first value displayed is the average value of all stored mea-

surements together with the WHO classification.

The total number of stored measurements is displayed at the same time, e.g. "A30". Repeated pressing of „M“ the individual measurement values are then displayed in sequence.

 This quick guide does not replace the detailed information on the use and safety of the blood pressure monitor featured in the user instructions.

Please be sure therefore to read the user instructions!

# Introduction

Dear Customer,

We are delighted that you have decided to purchase a boson blood pressure monitor. The boson brand is synonymous with optimal quality and precision and is also number 1 among professionals: 77% of all German doctors currently use boson blood pressure monitors. The experience gained from decades of use by professionals has also been incorporated into all the patient monitors for home measurement.

This instrument has also passed our strict quality control procedures and is your reliable partner for monitoring your blood pressure values.

Please read these user instructions before using the instrument for the first time, as correct blood pressure readings

can only be obtained if the instrument is operated correctly.

In these user instructions, „➡“ means an action by the user.

To get help concerning usage or maintenance please contact your dealer or the manufacturer (contact details can be found on the back cover of this manual).

An unexpected operating condition or an incidence which has worsened the health condition or could have worsened the health condition should be reported to the manufacturer immediately.

If the instrument is to be sold, ensure that these user instructions are enclosed.

## Introduction

If the monitor is used in medical practice (in accordance with the legislation governing the operation of medical devices), it should be calibrated at regular intervals (see page 72).

GB

The bosomedicus family 4 blood pressure monitor is suitable for patients with an upper arm circumference of between 22 and 48 cm. The monitor is not suitable for neonates.

Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 11 feet should be kept from such devices.

## Purpose

Non-invasive recording of systolic and diastolic blood pressure and pulse frequency in humans.

The instrument can be used with any arm circumference as shown on the corresponding cuff.

## Blood pressure values

To determine the blood pressure, **two values** need to be measured:

- **Systolic (upper) blood pressure:**

This is produced when the heart muscle contracts and blood is forced into the blood vessels.

- **Diastolic (lower) blood pressure:**

This is present when the heart muscle expands and fills with blood again.

Blood pressure measurement values are expressed in mmHg (mm of mercury).

The World Health Organization (WHO) has defined the following cut-off values for assessing blood pressure:

	systolic	diastolic
High blood pressure (hypertension)	over 140 mmHg	over 90 mmHg
Borderline normal blood pressure	130 to 139 mmHg	85 to 89 mmHg
Normal blood pressure	120 to 129 mmHg	80 to 84 mmHg
Optimum blood pressure	up to 119 mmHg	up to 79 mmHg

## Blood pressure values

---

Whether medical treatment is necessary depends not only on the blood pressure, but also on the patient's risk profile. Please contact your GP if one of the values (systole, diastole) consistently exceeds the cut-off point for high blood pressure.

GB

Self-measurement of blood pressure often produces somewhat lower values at home

than at the doctor's. The German Hypertension League has therefore defined lower cut-off values for home measurements:

Home measurement:	135/85 mmHg
Measurement in the doctor's surgery:	140/90 mmHg

## Starting up

The device must be installed and put into operation in accordance with the instructions in this manual.


Illustration 1

➔ Use only high quality leak-proof batteries of the correct specifications (see "Technical data" P. 69).

**⚠** **Never** mix old and new batteries or different makes.

**⚠** If the batteries are not inserted correctly, the monitor will not work and this may cause the batteries to heat up and leak, damaging the monitor.

### ➔ Inserting the batteries

The battery compartment is located on the underside of the monitor. Insert batteries in the battery compartment as shown in illustration 1.

➔ Remove the batteries if the monitor is not to be used for any length of time.

## Starting up

The monitor has a battery status display (see P. 40)

Battery fully charged.

Battery partially charged.

Replace battery soon.

Flashing! No further measurements possible.  
Replace battery.

### Care for the environment

Do not dispose of used batteries and charger units in domestic waste. You can take these to collection sites for used batteries or special waste. Contact your local authority for information.

## Set Time and Date

Your device bosomedicus family 4 has a clock function.

To select date and time setting, press the <SET> button. Modification of the parameters are done with the <START> - button.

➔ Press <SET> button for about 2 seconds, release - setting of date (year flashes), set year with START button.

➔ briefly press <SET> - setting of date (month flashes) , set month with START button

## Set Time and Date

➡ briefly press <SET> -  
setting of date (day flashes) , set day with  
START button

➡ briefly press <SET> -  
setting of time (hour flashes) , set hour  
with START button

➡ briefly press <SET> -  
setting of time ( minutes flashing) , set  
minute with START button

➡ briefly press <SET> -  
Setting of time mode (12h/24h/off flashes),  
set time mode with START button  
If you select 'off', no time will be displayed  
and the automatic switch-off period of the  
device will be extended to 60 minutes.  
N.B., this significantly increases battery

consumption, and so is recommended only  
when running the device on mains power.

➡ Briefly press <SET> , power is swit-  
ched off

## General instructions for self-measurement

1. Fluctuations in blood pressure are quite normal. Even when measurements are repeated, there can be marked differences. Single or irregular measurements do not provide reliable information about the actual blood pressure. A reliable evaluation is only possible if regular measurements are taken under comparable conditions and the measurement values entered in the blood pressure record card.
2.  Self-testing does not mean self-treatment. Do not on **any** account change the dose of medicine prescribed by your doctor of your own accord.
3. Irregular heartbeats can affect the accuracy of the monitor or result in incorrect measurements (see page 65).
4. Incorrect measurements can also occur in patients with a weak pulse (possibly in patients wearing a pacemaker, for example). The blood pressure monitor itself has no effect on the pacemaker.
5. If you are pregnant, consult your doctor before using the monitor.
6. Always take the measurement in a calm and relaxed state.  
You are recommended to check your blood pressure twice a day – in the morning on getting up and in the evening once you have relaxed after work.
7. Always measure your blood pressure from the arm with the higher blood pressure values (where not required otherwise (see page 54)).

## General instructions for self-measurement

---

To do this, measure the blood pressure first in both arms and thereafter always from the arm with the higher blood pressure.

## Preparing to measure your blood pressure

### Note on cuff sizes:

 The monitor may only be used with the following cuff types. These should be chosen according to the upper arm circumference printed on them.

GB

Type	Arm circumference	Order number
CA01	22 – 32 cm	143-4-764
CA02	32 – 48 cm	143-4-757
CA04	22 – 42 cm	143-4-765*

\* Universal cuff supplied as standard.

Attach the cuff to the monitor by pressing the plug (click) into the socket on the left side of the unit. (See illustration 2).


Illustration 2

## Preparing to measure your blood pressure

Observe the following instructions when measuring your blood pressure:

Avoid smoking or drinking coffee for one hour before measuring your blood pressure.

Sit comfortably for the blood pressure measurement. Support your back and arms. Do not cross your legs. Place your feet flat on the floor.

Relax for 5 minutes before taking a measurement.

Do not move during the measurement.

### Attaching the cuff

The measurement should be taken from the bare upper arm.

In the case of tightly fitting outer clothing, ensure that the blood supply to the arm is not cut off when the sleeve is rolled up (if necessary remove the article of clothing before taking the measurement).

➔ Open out the cuff into the shape of a ring and slide it over the upper arm until the bottom edge of the cuff is about 2-3 cm above the elbow. The cuff must be positioned so that the marking is over the artery (see illustration 3).

## Preparing to measure your blood pressure


Illustration 3

GB

⚠ Do not attach the cuff over wounds as this can cause further injuries.

⚠ Ensure that the cuff is not attached to an arm where arteries or veins are undergoing or have undergone medical treatment (e.g. shunt).

⚠ For women who have had a breast amputated, do not attach the cuff to the arm on the amputated side of the body.

⚠ During the measurement, medical apparatus that are being used on the same arm at the same time, may malfunction.

⚠ The cuff must not be too tight. There should be enough space to fit roughly two fingers between arm and cuff.

⚠ The metal ring should **never** lie over the artery (see illustration 3) as this may produce a false reading.

➔ Grasp the end of the cuff that feeds through the metal ring and wrap it around the outside of your arm. Attach the Velcro fastener by applying gentle pressure.

## Preparing to measure your blood pressure

- ➔ Place your arm with the cuff attached on the table, holding it relaxed and slightly bent so that the cuff is level with your heart.


Illustration 4

Do not talk during the measurement.

**⚠** Check that the air tube is not kinked during the measurement. This could cause congestion of the blood flow and subsequent injury.

**⚠** The blood pressure measurement should not prevent the flow of blood for an unnecessary length of time (> 2 minutes). If the monitor fails to function correctly, remove the cuff from the arm.

**⚠** Excessive frequent measurements can cause injury by restricting blood flow.

## Select an user

Your device boso-medicus family 4 has 60 memories for 4 persons.

Pressing the <SET> button with the power switched off allows to select one of the users.

GB

Then press <START> button to determine the selected person.

## Measuring your blood pressure

➔ Start the measurement with the start button

⚠ Now hold the arm absolutely still and do not talk.

All the items on the LCD display appear briefly to indicate that the monitor is now ready for use (see illustration 5).


Illustration 5

## Measuring your blood pressure

The pump starts to inflate the cuff.

The monitor has an intelligent automatic inflation system for gentle inflation to the correct cuff pressure.

The increasing cuff pressure is displayed.

When the **required pressure** is reached, the pump switches off and the air is released slowly from the cuff.

### **Interrupting a measurement:**

A measurement can be interrupted at any stage by pressing the start button; the cuff deflates automatically.

During the measurement phase the descending cuff pressure and the flashing  symbol are displayed (see illustration 6).

Please keep your arm still if the symbol  appears while a measurement is being taken.



Illustration 6

## Measurement display

After the end of the measurement, the in-built valve opens automatically and the cuff deflates rapidly.

The measured blood pressure values (systole, diastole, pulse) and time and date are displayed. see illustration 7).

Alternating time and date are displayed.

Illustration 7

The measurement is stored automatically (except for guest mode, see below). If the memory store is full (60 measurements), the **oldest** measurement is **deleted** and the current one stored.

### Note on guest mode:

If a measurement value is not be kept, press „M“ while the measured values are still on display.

Invalid measurements with an “Err” message (see page 64) are not stored.

## Measurement display

If the symbol "♥" appears after the measurement, it is recommended that the measurement is repeated, holding the arm completely still. If the "♥" symbol appears again even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

The "👤" symbol indicates that the arm moved during measurement. You are advised to take another measurement.

Blood pressure is a dynamic parameter and can be affected by the patient's position, e.g. sitting, standing, lying, moving, before or during the measurement, physical condition (stress, disease, etc.).

- ➔ If the measurement result is obviously wrong, repeat the measurement.

The unit switches itself off automatically after about 1 min. To repeat a measurement, press the start button again **briefly**.

Allow at least two minutes between measurements.

- ➔ If no further measurement is required, remove the cuff from the arm.

# Completing the blood pressure record card

## Retrieving the average value and WHO classification

➔ To recall the stored press the M key with the device switched off.

The average values of all stored measurements (systolic, diastolic and pulse) are displayed.

GB

➔ Pressing the <SET> button during the memory-recall allows to select one of the users.

If there is no measurement in the memory, "0" is displayed for systole, diastole and pulse.

If measurements were recorded with irregular heart beat, there is an additional symbol showing the prevalence of cardiac arrhythmias:

appears when up to 25% of the measurements have irregular heart beat.

appears when between 25% and 50% of the measurements have irregular heart beat.

appears when more than 50% of the measurements have irregular heart beat.

# Completing the blood pressure record card

## Retrieving the average value and WHO classification

A dark rectangle appears on the left of the screen (see illustration 8), indicating the level of the average blood pressure value. Use the accompanying evaluation scale (see illustration 9) to identify the range within which your blood pressure falls under the WHO classification (see page 45).

➔ For your information, attach the sticker (WHO classification) supplied, to the back of your unit.


Illustration 8

### WHO classification

- high blood pressure (hypertension)
- borderline normal blood pressure
- normal blood pressure
- optimum blood pressure

Illustration 9

# Displaying individual stored values and clearing the memory store


➔ Recall the last measurement stored in the memory by pressing „M“ again after having recalled the average value. The number of the measurement appears in the display (see illustration 10).

➔ Repeated pressing of „M“ will now successively display all measurements contained in the memory as described above.

The corresponding measurement value appears after 3 seconds (see illustration 11).

➔ Pressing the START button switches off the device.

GB

Systole  
mmHg

Diastole  
mmHg

Pulses/min.

Alternating time and date  
are displayed.

If no button is pressed, the device automatically switches off after 60 seconds.

## Clearing the memory store

Press „M“ and hold for 5 seconds until the symbol  starts to flash. The memory of the currently selected user is erased.

Illustration 10

Illustration 11

## Mains operation

---

The mains connection socket is located at the rear of the monitor.


Where required, use only the bosco power supply unit (**order number 410-7-150**). This power supply unit produces a rectified output of the correct polarity. Other commercially available power supply units can cause damage to the electronic components, which will invalidate the warranty. There is also a risk of fire with these other power supply units.

## Error messages

If any problems occur during the measurement, the display will show an error message instead of a blood pressure reading.

Meaning of error messages:

	Cause of <b>Err</b> message:	Troubleshooting:
 	Excessive variation in cuff pressure during measurement	Hold arm still.
	No valid pulse readings	Check position of cuff (see page 53) and repeat measurement.
	Systole – diastole difference too small: Systole – diastole $\leq 10$ mmHg	Check position of cuff (see page 53) and repeat measurement.
	Incorrectly inflated	Cuff possibly too loose, attach cuff more firmly (see page 53)

## Error messages

Cause of <b>Err</b> message:		Troubleshooting:
	The pulse value could not be determined correctly.	Check position of cuff and repeat measurement.
	Internal Error	Remove the batteries and briefly press START button. Put the batteries back into the monitor. If the error persists, return the instrument for repairs to the manufacturer.
		
	"  Display:	Troubleshooting: Repeat measurement recommended, holding the arm absolutely still. If the "  symbol appears again, even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

## Cleaning and Disinfection

---

### **Cleaning Monitor:**

Use a soft, dry cloth to clean your monitor.

### **Cleaning Cuff:**

Small stains on the cuff can be removed carefully with a proprietary washing-up liquid.

The logo consists of the letters 'GB' in a bold, sans-serif font, enclosed within a white circle. This circle is centered on a dark grey rectangular background.

### **Disinfection Cuff:**

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant antiseptic liquid (Schülke & Mayr). To disinfect the velcro of the cuff, we recommend spray disinfection.

In particular, if the device is used on multiple users is to pay attention to regular cleaning and disinfections of the cuff.

## Warranty Conditions

---

We give 3 years warranty from the date of purchase. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole unit, but only to the replaced components.

Excluded from the warranty are parts subject to normal wear and tear (e.g. cuff), transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty.

No claims for damages against us are substantiated by the warranty.

In the case of justified warranty claims the device has to be sent along with the original invoice to:

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## Customer service

### Customer service:

Warranty and repair work must be undertaken by trained and authorised personnel. Do not modify this equipment without authorization of the manufacturer.

Send the monitor, carefully packaged and with sufficient postage, to your authorised dealer or directly to:

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## Disposal

Do not dispose of monitors and batteries in the domestic waste.

At the end of its lifespan, the monitor must be taken to a collection site for obsolete electronic items.

Please note that batteries and rechargeable batteries must be disposed of separately (local authority collection site).

Expected operational lifetime of the equipment: 5 years (using the unit 6 times a day).  
Expected operational lifetime of the cuff: 2 years (using the unit 6 times a day).

## Technical data

---

<b>Measurement principle:</b>	oscillometric
<b>Measurement range:</b>	Systolic: 60 to 279 mmHg Diastolic: 40 to 200 mmHg Pulse: 40 to 180 per minute
<b>Cuff pressure:</b>	0 to 299 mmHg
<b>Memory store:</b>	4 x 60 measurements
<b>Display:</b>	LCD
<b>Operating conditions:</b>	environmental temperature +10°C to +40°C Relative humidity 15 to 85%
<b>Storage conditions:</b>	environmental temperature -20°C to +60°C Relative humidity 10 to 95%
<b>Power supply:</b>	DC 6 V (4 x 1.5 V Mignon IEC LR 6 alkaline manganese batteries) Alternative special option: DC 6 V power supply unit, Order No. 410-7-150

## Technical data

---

<b>Typical battery life:</b>	700 measurement cycles (depending on inflation pressure and frequency of use).
<b>Battery check:</b>	symbol display in LCD window
<b>Weight:</b>	282 g without batteries
 <b>Dimensions (WxHxD):</b>	140 mm x 60 mm x 106 mm
<b>Classification:</b>	Protection class II (□)
<b>Applied part:</b>	Type BF (⚡)
<b>Protection against solid objects:</b>	IP20
<b>Clinical test (DIN 58130):</b>	Accuracy complies with the requirements of EN 1060 Part 3

## Technical data

---

**Maximum deviation of cuff pressure measurement:**  $\pm 3$  mmHg

**Maximum deviation of pulse rate display:**  $\pm 5$  %

**Applicable standards:** EN 1060, Part 1: "Non-invasive sphygmomanometers – General requirements" and Part 3: "Supplementary requirements for electro-mechanical blood pressure measuring systems".

# Calibration checks - Testing instructions

## (for trained specialist staff only)

Calibration checks – every 2 years at the latest – must be performed by the following institutions or persons:

- Manufacturer
- Trade measurements authorities
- Persons who fulfil the requirements of the legislation governing the operation of medical devices.

### A) Function testing

Function testing of the monitor can only be carried out on a person or with a suitable simulator.

### B) Testing of pressure circuit integrity and deviation of pressure display

#### **Note:**

- a) When the pressure in measurement mode increases above 299 mmHg, the quick release valve is triggered and the pressure circuit opens. When the pressure in test mode increases above 320 mmHg, the measured value display flashes.
- b) For the calibration checks a special test plug must be used. The test plug is available for a fee from the manufacturer.  
A pump ball must also be introduced into the pressure circuit.

# Calibration checks - Testing instructions

## (for trained specialist staff only)

### Testing

- 1.) Remove batteries.
- 2.) Hold start button pressed down and insert batteries.
- 3.) Release start button.  
The monitor is now ready for testing; the current pressure is displayed in the SYS and DIA fields.
- 4.) Perform test for deviation of pressure display and pressure circuit integrity (observe setting time for cuff – at least 30 seconds) in the usual manner.

5.) Switch off monitor by pressing the START button again.

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**BOSCH + SOHN GmbH u. Co. KG**